

Recommended protocol for CD8

Obtained in run 72

06 Sep 2024

Immunostainer

Type: intelliPATH

Primary antibody

Clone: C8/1779R
Producer: Biocare Medical
Product no. / lot no.: API3219 / TBD
Diluent: None
Diluent producer / no.: (None) / None
Dilution factor: 1:0
Incubation time / temperature: 30 min. / 25°C

Epitope retrieval, HIER

Device: Pressure Cooker
Buffer: Biocare Medical Diva Decloaker
Heating time at max. temp.: 15 min.
Maximum heating temp.: 110°C

Visualization system

Producer: Biocare Medical
Product / no: MACH 4 HRP-Polymer / MRH534
Linker: None
Incubation time polymer: 30 min.
Incubation temperature: 25°C

Chromogen

Producer: Biocare Medical
Product / no: intelliPATH DAB / IPK5010
Incubation time / temperature: 5 min. / 25°C
Enhancement: None

Disclaimer:

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.