

## Recommended protocol for OCT3/4

Obtained in run 53

13 Feb 2018

### Immunostainer

Type: Manual

### Primary antibody

Clone: NRG1.1  
Producer: Diagnostic BioSystems  
Product no. / lot no.: Mob 464-05 / W053  
Diluent: Diamond - Antibody Diluent  
Diluent producer / no.: Cell Marque / 938B  
Dilution factor: 1:50  
Incubation time / temperature: 900 min. / 4°C

### Epitope retrieval, HIER

Device: PT-link / PT-module  
Buffer: Dako TRS High pH (3-1)  
Heating time at max. temp.: 20 min.  
Maximum heating temp.: 98°C

### Visualization system

Producer: Dako  
Product / no: EnVision FLEX+ / K8012/DM822  
Linker: Linker, Mouse  
Incubation time linker: 10 min.  
Incubation time polymer: 10 min.  
Incubation temperature: 37°C

### Chromogen

Producer: Dako  
Product / no: DAB+ / DM827 / DM823 / K800  
Incubation time / temperature: 5 min. / 37°C  
Enhancement: None

### Disclaimer:

NordiQC makes every attempt to provide accurate and up-to-date information, yet NordiQC does not make any claim or warranty regarding the accuracy of the provided information nor does it represent that the contents of the web site and protocols reflect the most recent developments in immunohistochemistry at any point in time.